







**6. Serum Electrophoresis**

	Test	Method	Specification	Result
6.1	Albumin	Agarose Gel Plate	As reported	<b>18 g/L</b>
6.2	Alpha 1 Globulins	Agarose Gel Plate	As reported	<b>5 g/L</b>
6.3	Alpha 2 Globulins	Agarose Gel Plate	As reported	<b>8 g/L</b>
6.4	Beta Globulins	Agarose Gel Plate	As reported	<b>4 g/L</b>
6.5	Gamma Globulins	Agarose Gel Plate	As reported	<b>1 g/L</b>
6.6	Electrophoresis Pattern	Agarose Gel Plate	Typical	<b>Typical</b>

**7. Antibody testing**

	Test	Method	Specification	Result
7.1	Akabane Virus	ELISA	Not detected	<b>Not detected</b>
7.2	Bluetongue Virus	ELISA	Not detected	<b>Not detected</b>
7.3	Bovine Herpesvirus 1	VNT	Not detected	<b>Not detected</b>
7.4	Bovine Viral Diarrhoea	VNT	As detected	<b>Titre 32</b>
7.5	Enzootic Bovine Leukosis	ELISA	Not detected	<b>Not detected</b>
7.6	Parainfluenza 3	ELISA	Not detected	<b>Not detected</b>

Measles, mumps, and rubella is not infect *bos taurus*, and rabies virus is not present in Australia.

The products described above were derived from animals in Australia and the products were treated solely within Australia. Foot-and-mouth disease and rinderpest do not exist in Australia.

The serum is intended for research or pharmaceutical/technical purposes only. The country of import/origin is free from TSE group of diseases and the product is free from pathogens especially Bovine Viral Diarrhoea (BVD).

**Signed** \_\_\_\_\_  
 QC